

# Eliason Lab Respo Cleaning Guidelines

April 10 2020

K. Kraskura, E. Hardison, T. Dressler, J. Van Wert


## **Data Analysis Pipeline:**

- 1) Visualize Data
  - a. Plot all MO2 data using Krista's code
  - b. Visually inspect all MO2 plots (all files - background, MMR, SMR)
- 2) Clean Data
  - a. Identify all clean MO2 data with linear decline - keep
  - b. Identify any plots with minor issues - correct
  - c. Identify erroneous MO2 plots - remove
- 3) Background
  - a. Assess if background is present
  - b. If background is present, correct all MO2 slopes for background respiration (see Rosewarne et al 2016 for guidance)
  - c. We generally assume a linear increase between the "before" and "after" background and apply that to the slopes.
- 4) Plot your data and Calculate values
  - a. ALWAYS plot MO2 vs time and temperature vs time. ESSENTIAL for every fish.
  - b. Calculate SMR, RMR, MMR, EPOC, time to MMR50, MMR 75
  - c. Many ways to do this – talk to Erika and each other (i.e. which sliding measurement window; which SMR value; which EPOC threshold)
- 5) Assess for scaling
  - a. Plot your data – SMR vs body mass; MMR vs body mass
  - b. Adjust for scaling if needed (see above).

## Notes before Cleaning:


### 1. Check for leaks

- ❖ Plot flush plot ON
- ❖ Are L shapes consistent across files? If yes, you might have a leak.




1. Note: Do not cut your measurement down to below your “sliding window” length.
2. For cleaning:  $R^2$  below 0.7 are immediately rejected-- no cleaning
3. If you have any measurements where you have a note that you were messing with the chamber (relieving a bubble, etc) or where data is missing from the measurement (low signal in probe can cause this), exclude the measurement

# Decision Tree - Oxygen Visualization




# Scatter





if there are linear decreases with a slightly reduced  $R^2$  due to scatter, this is likely because the measurement cycle was too short or because the O<sub>2</sub> probe is old and signal is poor. If these look like clean linear oxygen consumption apart from the points being spread more than usual, adjust your  $r^2$  threshold accordingly (nothing below 0.7) for all your files.

# Z-Shape - three part change




1. IS THIS SETUP RELATED? (e.g. huuge fish in small chamber, where every data point is super influenced by actual consumption of the fish) - if yes, don't clean trust the avg regression.
2. If this is an anomaly— we propose excluding this measurement. Could be due to a mixing or pump problem. E.g. Recirculation pump has seized and water is not flowing consistently past O<sub>2</sub> probe. Or there is a bubble in T joint adjacent to probe. Data cannot be trusted

**L- shape** - Bottom of L must be  $>1/4$  of the trace to be classified as an L



Cut the bottom of the L off




# Tiny Tails - ends deviate from linearity (small $\tau$ or small $L$ )



Chop off end



**7's** - top of measurement is flat and  $>1/4$  of the file




Cut the top of the 7 (reverse L)

# Double Slope - two clearly defined decreasing slopes

You have a choice on how to handle this. **You just need to be consistent** in YOUR analysis and it needs to make sense for YOUR goals of your experiment. So, make a decision and stick with it for your entire study.





Options:

- 1) If it is rare, and you have lots of data, just delete this file.
- 2) You can take the longer slope (rationale = most representative of this measurement period), if it meets R2 cut-off, then use it.
- 3) You can take the shallower slope (rationale = we are trying to measure minimum metabolism estimates and this is real data representative of the fish), if it meets the R2 cut-off.
- 4) (Less likely, but possible) If you are trying to estimate MMR, see just how fast the steep slope is in case it meets your MMR criteria. If it does, use it for your MMR.



# Concave and convex (Note that this is for individuals that do not have a leak in the respo)

\*Common during EPOC




**Common procedure:** Convex traces can be real, fish is slowly increasing its O<sub>2</sub> consumption. I assign an R<sup>2</sup> cut-off for this single measurement (recommend 0.05 higher than your global cut-off) and keep it if it's within that threshold because it is probably good data.

**NOTE:** Evaluate how frequent these traces are. Could make a personal judgment call for an alternate options:

1. For long measurements (trace 1): chop that data into smaller time periods. Min 5-8 = slope 1, Min 8-11 = slope 2, Min 11-14 = slope 3. You have a recovering fish that has very different metabolic rates over this 10 minute period and you are losing really interesting discrimination power by averaging over this entire time period.
2. If it's rare and not at critical time point, delete.



Either handling choice is ok, but must stay consistent across all trace analysis for any given study.

# Zig Zags, blips, and blops in the middle



Could be caused by temperature probe; let the R2 filter

# Step down (lightning bolt)





Look at temp plot. Is the plateau due to temp?



Yes → Trim out the plateau

No → Reject trace

# Specifically for MMR

1. Visualize your data and look for blips in your MMR measurement
2. For files where there ARE blips, double check where the sliding window is calculating MMR.
3. If the sliding window with the steepest slope is over that blip, remove the blip from your analysis and re-run your MMR file.



Problem	Potential issue	Protocol	Example
Sudden, rapid decline in [O <sub>2</sub> ]	poor mixing or issue with recirc pump. If the fish truly respired that much and decreased the DO levels in the chamber that low, the flush pump would not be able to return the DO back up to full saturation in one single flush cycle	Remove the slope(s)	
Consistent non-linearity	Examine all trials. Cross-check with temp. Respirometer might have a leak and data is unreliable.	X	
Overshoot	Sudden large drop in [O <sub>2</sub> ] and flush then overcompensates. Common in animals with high metabolic rates (i.e. salmon)	Remove the slope(s)	